

## Inverted spinning disk confocal



**Fast confocal** - Equipped with a Yokogawa CSU-X1 spinning disk unit for the ultimate in high speed and low phototoxicity for live cell imaging.

Microscope	Nikon Eclipse Ti Inverted
Imaging System	Yokogawa CSU-X1 disk head and Andor Neo sCMOS camera
Objectives	4x, 10x, 20x (Plan Apo VC), 40x (ELWD), 60x oil and 100x oil.
Laser wavelengths	405nm, 488nm, 561nm, 640nm